

Analog Dry Bath Incubators



Analog Dry Bath Incubator One Block
Model 112001



Analog Dry Bath Incubator Two Block
Model 112002



Analog Dry Bath Incubator Four Block
Model 112004

- One, two and four block capacity models
- Amber Pilot Light
- Hydraulic Thermostat
- Dynamic temperature range
- Excellent value
- Blocks are interchangeable with Boekel Digital Models

Supplied with Block Extractor

Analog Dry Baths Specifications			
One Block Capacity		Model 112001	Model 112001-2
Dimensions	in	8.0 in x 9.8 in x 3.5 in	
W x L x H	cm	20.3 cm x 24.9 cm x 8.9 cm	
Net Weight		4.0 lbs	
Shipping Weight		4.5 lbs (2.0 kg)	
Temperature Range		up to 150°C	
Uniformity		+/- 0.50°C @ 37°C	
Stability		+/- 0.45°C @ 37°C	
Electrical		115 VAC, 50/60 Hz, 100 W	230 VAC, 50/60 Hz, 100 W (CE)
Two Block Capacity		Model 112002	Model 112002-2
Dimensions	in	8.0 in x 9.8 in x 3.5 in	
W x L x H	cm	20.3 cm x 24.9 cm x 8.9 cm	
Net Weight		4.5 lbs	
Shipping Weight		5.0 lbs (2.0 kg)	
Temperature Range		up to 130°C	
Uniformity		+/- 0.50°C @ 37°C	
Stability		+/- 0.45°C @ 37°C	
Electrical		115 VAC, 50/60 Hz, 150 W	230 VAC, 50/60 Hz, 150 W (CE)
Four Block Capacity		Model 112004	Model 112004-2
Dimensions	in	8.0 in x 9.8 in x 3.5 in	
W x L x H	cm	20.3 cm x 24.9 cm x 8.9 cm	
Net Weight		6.0 lbs	
Shipping Weight		6.5 lbs (2.5 kg)	
Temperature Range		up to 130°C	
Uniformity		+/- 0.50°C @ 37°C	
Stability		+/- 0.50°C @ 37°C	
Electrical		115 VAC, 50/60 Hz, 300 W	230 VAC, 50/60 Hz, 300 W (CE)