


## Heating Block Modules

Boekel offers the complete spectrum of Block Modules. These precision-machined solid aluminum blocks can accommodate microcentrifuge tubes, microplates and test tubes up to 25 mm in diameter. We have added blocks for use with 96 well 0.2 ml conical microtubes, 1.5 ml "Dolphin" tubes used for extractions, 30 ml nickel crucibles. Each block has a well for thermometer (not included) and a block extractor hole.

Block Modules are interchangeable between analog and digital dry bath incubator models. All blocks can also be used in the Boekel MiniFridge II for cooling applications.



These anodized blocks resist corrosion and feature alpha/numerical labeled rows and columns for easy vessel identification.

**Microplate Block Modules** – Applications include incubations for EIA's, ELISA's and monoclonal antibody

screening. The #110051 holds one round bottom microplate. For microplate cooling applications, The Microplate Block can be used in the Boekel MiniFridge II benchtop cooler. The #110096 96 well microtube block is designed to handle plates, tube strips or 0.2 ml tubes. Microplate Blocks fit both Digital and Analog Dry Baths.

### Test Tube Block Modules

	110006	110010	110013	110016	110020	110025	110035	110040	110045	110048
Dimensions	in				3 x 3.75 x 2					
W x L x H	cm				7.6 x 9.5 x 5.1					
Net Weight	1 lb (0.45 kg)									
Shipping Weight	2.25 lb (1.0 kg)									
Capacity	35 - 6 mm tubes	20 - 10 mm tubes	20 - 13 mm tubes	12 - 16 mm tubes	6 - 20 mm tubes	6 - 25 mm tubes	3 - 25 mm tubes; 5 - 13 mm tubes; 6 - 6 mm tubes	24 - 1.5 mL tubes	24 - 0.5 mL tubes	48 - 0.2 mL tubes

### Miscellaneous Block Modules

	110051	110096	110099	
Dimensions:	in	3.75 x 5.25 x 1.75	3.75 x 6 x 1	3 x 3.75 x 2
W x L x H	cm	9.5 x 14.3 x 4.4	9.6 x 15.3 x 2.5	7.6 x 9.5 x 5.1
Net Weight		2.25 lb (1.1 kg)	2.25 lb (1.1 kg)	1.25 lb (0.5 kg)
Shipping Weight		3 lb (1.4 kg)	3 lb (1.4 kg)	2 lb (0.9 kg)
Capacity		1 - 96 well microplate	1 - 96 well 0.2 mL microplate; 0.2 mL strips or 0.2 mL tubes	24 - 1.5 mL "Dolphin" tubes

### Options

Acrylic Plastic Cover for one and two block dry baths  
Acrylic Plastic Cover for four block dry baths

### Model

110002LID  
110004LID



855 Pennsylvania Boulevard • Feasterville, Pennsylvania, USA 19053  
Phone: 215-396-8200 or 800-336-6929 • Fax: 215-396-8264  
Email: boekel-info@boekelsci.com • Web site: <http://www.boekelsci.com>