

BOEKEL SLIDE STORAGE CABINET

Model 141000

FEATURES:

- Six cabinet drawers hold up to 5000 standard 1 inch x 3 inch microscope slides
- Cabinets can be interlocked up to 10 units high
- Plastic drawer dividers can be removed to store objects up to 5 cm (2 in.) wide
- Drawers slide smoothly and easily
- Back of drawer acts as handle for transportation
- Durable cold rolled steel with baked on epoxy finish

SUPPLIED WITH:


- Sponge block spacers
- Identification labels

ACCESSORIES:

- Cabinet Base with Rubber Feet
- Drawer Only

MODEL

141001
141002



Shown with Two Microscope Slide Storage Cabinets (Model 141000) and One Cabinet Base Model (141001)

Slide Storage Cabinet and Cabinet Base Specifications

	Slide Storage Cabinet Model 141000	Cabinet Base Model 141001
Overall Dimensions: (Width x Depth x Height)	40 cm x 48.3 cm x 12.7 cm (15-3/4 in x 19 in x 5 in)	40 cm x 47.6 cm x 7 cm (15-3/4 in x 18-3/4 in x 2-3/4 in)
Net Weight	7.7 kg (17 lb)	1.3 kg (2-3/4 lb)
Shipping Weight	9.1 kg (20 lb)	1.8 kg (4 lb)

NOTE: Cabinets are stacked by built-in tabs in the bottom of each unit. Because of this, cabinets will "rock" and be unstable if placed on a bench-top without the use of Cabinet Base 141001 underneath the first cabinet. Subsequent cabinets may then be stacked on top of the first cabinet. See photo.

BOEKEL MICROSCOPE SLIDE DISPENSER

Models 1415

FEATURES:

- Stainless steel slide dispenser holds up to 54-3 in. x 1 in. slides
- Adjustable to accommodate slides of different thickness
- Rubber cushion on bottom prevents unit from slipping and scratching the benchtop



Boekel Microscope Slide Dispenser Specifications

Overall Dimensions: (Width x Length x Height)	7.6 cm x 12.7 cm x 10.2 cm (3 in x 5 in x 4 in)
Net Weight	0.5 kg (1 lb)
Shipping Weight	0.9 kg (2 lb)

BOEKEL MICROSCOPE SLIDE DISPENSER has a 2 year warranty on parts and workmanship!

Boekel Scientific manufactures life science and general lab products, including hybridization equipment, benchtop coolers and freezers, platform mixers, microplate equipment, incubators, dry bath incubators, clinical products, water baths and circulators, ultrasonic cleaners, aspirator pumps, ovens, desiccators, pipet handling equipment and lab ware..



855 Pennsylvania Boulevard • Feasterville, Pennsylvania, USA 19053
Phone: 215-396-8200 or 800-336-6929 • Fax: 215-396-8264
Email: boekel-info@boekelsci.com • Web site: <http://www.boekelsci.com>