

Boekel Slide Storage Cabinets

Model 141000

Features:

- Six cabinet drawers hold up to 5000 standard 1 inch x 3 inch microscope slides
- Cabinets can be interlocked up to 10 units high
- Plastic drawer dividers can be removed to store objects up to 5 cm (2 in) wide
- Drawers slide smoothly and easily
- Back of drawer acts as handle for transportation
- Durable cold rolled steel with baked on epoxy finish



Supplied With:

- Sponge block spacers
- Identification labels

Accessories:

- Cabinet Base with Rubber Feet, Part Number: 141001
- Drawer Only, Part Number: 141002

	Slide Storage Cabinet Model 141000	Cabinet Base Model 141001
Overall Dimensions (Width x Length x Height)	40 cm x 48.3 cm x 12.7 cm (15 3/4 in x 19 in x 5 in)	40 cm x 47.6 cm x 7 cm (15 3/4 in x 18 3/4 in x 2 3/4 in)
Net Weight	7.7 kg (17 lb)	1.3 kg (2 3/4 lb)
Shipping Weight	9.1 kg (20 lb)	1.8 kg (4 lb)

Note: Cabinets are stacked by utilizing built-in tabs in the bottom of each unit. Because of this, cabinets will "rock" and be unstable if placed on a bench-top without the use of a Cabinet Base, Part Number: 141001, underneath the first cabinet. See photo.

Boekel Wax Dispenser

Model 145600

Features:

- 5 Liter capacity, Approximately 1.25 Gallon
- Heated valve ensures a consistent flow of molten wax
- Stainless steel outer casing
- Aluminum inner reservoir provides excellent thermal transfer and temperature control
- Dispenser valve allows precise control of flowing wax and can be locked in the open position if desired
- Inner reservoir screen helps to prevent large particulate from clogging the dispenser valve
- Rapid warm-up with an advanced safety over-temperature cut-off

Wax Dispenser	
Model Number	145600 (115 VAC), 145600-2 (230 VAC) CE
Overall Dimensions	21.6 cm x 40.6 cm (8 1/2 in x 16 in)
Temperature Range	40° C to 75° C +/- 2° C
Approvals	Certified to UL61010-1 and CSA C22.2
Electrical	115 VAC 50/60 Hz 250 W



855 Pennsylvania Boulevard • Feasterville, PA 19053 USA
 800-336-6929 • Phone: 215-396-8200 • Fax: 215-396-8264
 E-mail: boekel-info@boekelsci.com • Web site: www.boekelsci.com