

Boekel Lighted Tissue Flotation Baths

Models 145951 and 145702

Features:

- Eliminate wrinkles and distortion during preparation of paraffin embedded tissue sections
- L.E.D. illuminated black background with transparent PYREX® dish; allows specimen to be easily seen
- Digital microprocessor control with advanced safety over temperature protection
- Deluxe model includes Slide Dryer Block to assist in drying of slides and Tissue Orienter Block to smooth wrinkled sections

Supplied With:

- Removable PYREX® glass dish, Part Number: 904-0013

Accessories

- PYREX® glass dish, Part Number: 904-0013



145951 Deluxe Lighted Tissue Flotation Bath



145702 Standard Lighted Tissue Flotation Bath



Boekel Slide Dispenser Model 1415:

- Stainless steel slide dispenser holds up to 54 3 in. x 1 in. slides
- Adjustable thickness
- Rubber cushion on bottom
- More specs on www.boekelsci.com

	Standard Lighted Tissue Flotation Bath	Deluxe Lighted Tissue Flotation Bath
Model Number	145702 (115 VAC) 145702-2 (230 VAC)	145951 (115 VAC) 145951-2 (230 VAC)
Overall Dimensions (Width x Length x Height)	36.2 cm x 36.2 cm x 11.4 cm (14.25 in x 14.25 in x 4.5 in)	36.2 cm x 36.2 cm x 11.4 cm (14.25 in x 14.25 in x 4.5 in)
Dryer Block	NA	5.1 cm x 12.7 cm (2 in x 5 in)
Orienter Block	NA	5.1 cm x 5.1 cm (2 in x 2 in)
Net Weight		5.45 kg (12 lb)
Shipping Weight		6.8 kg (15 lb)
Temperature Set Range		Ambient +5° C to 65° C
Accuracy		+/- 1° C
Uniformity at 55° C		+/- 1° C
Dryer Block Temperature	NA	40° C +/- 3° C preset
Orienter Block Temperature	NA	55° C +/- 3° C fixed
Approvals		Certified to UL61010-1 and CSA C22.2
Electrical	145701 (115 VAC 50/60 Hz 272 W) 145702-2 (230 VAC 50/60 Hz 272 W) CE	145951 (115 VAC 50/60 Hz 272 W) 145951-2 (230 VAC 50/60 Hz 272 W) CE

Boekel Tissue Baths have a 2 year warranty on parts and workmanship