

Shake 'N Bake Digital Incubator

The Shake 'N Bake Incubator shelves provide a gentle up and down motion which is ideal for hybridization and wash procedures as well as conventional incubations. This unit features a digital temperature setting and display plus a PID temperature controller for rapid heat-up and excellent temperature stability.

Ordering Information

Part No.	Size	Electrical	Door Style
136400	0.7 cu. ft.	115 volt/50-60 Hz	See-thru
136400-2	0.7 cu. ft.	230 volt/50-60 Hz	See-thru

Model **136400 Shake N' Bake**

Chamber Volume	0.7 cu. ft.
Chamber Dimensions (in.) W x D x H (cm)	11-1/2 x 11-1/2 x 10-1/2 29.2 x 29.2 x 26.7
Temperature Range	Ambient +5°C to 70°C
Stability @ 37°C	±0.5°C at 65°C
Uniformity @ 37°C	±0.5°C
Shelves Included	2
Max No. Shelves*	2
Shelf Area	1.84 ft.2
Rocking Speed	8 to 25 strokes per minute
Max. Shelf Load	1 kg (2-1/4 lb.)
Dimensions (in.) W x D x H (cm)	12-1/2 x 14-1/4 x 16 31.8 x 36.2 x 40.6
Shipping Weight	28 lbs./12.7 kg



*For accessory shelf information contact Boekel at 800-336-6929 or visit our web site www.boekelsci.com/incubators



855 Pennsylvania Boulevard • Feasterville, PA 19053 USA
800-336-6929 • Phone: 215-396-8200 • Fax: 215-396-8264
E-mail: boekel-info@boekelsci.com • Web site: www.boekelsci.com

©Copyright, Boekel Scientific, 2009

Printed in USA

1109