



Boekel Tube Rotator

FEATURES

- Digital microprocessor control
- Adjustable speed provides gentle or vigorous mixing action
- Adjustable angled head offers mixing angles from 90 to 180 degrees
- Five carousel options provide greater versatility for the most tube mixing needs
- Timer features offers time-down or time-up counting feature
- Includes 260653 - 14 Place Tube Carousel for 50mL
- Power pac is UL, CSA and CE approved
- CE marked
- 2 year warranty

ACCESSORIES

260651	46 place tube carousel for 1.5/ 2.0 ml tubes
260652	28 place tube carousel for 15 mL tubes
260653	14 place tube carousel for 50 mL tubes
260654	63 place tube carousel for 18.5 mm cul. tubes
260655	63 place tube carousel for 14 mm cul. tubes



Boekel Orbital Shaker

FEATURES

- Compact with low profile and small footprint
- Smooth horizontal orbital motion for mixing in bottles, flasks, and beakers
- Soft start and noise-free operation
- Can be used in an incubator 4°C to 40°C
- Platform not supplied. See interchangeable platforms listed below.
- Power pac is UL, CSA and CE approved
- CE marked
- 2 year warranty

FLASK DANCER PLATFORMS

Catalog No.	Interchangeable Platforms
270290	6 x 250 mL flasks/400 mL beakers
270291	12 x 100 mL flasks/150 mL beakers
270292	Universal platform with adjustable bars
270293	Flat platform with rubber mat



Orbitron III™ Platform Rotator

FEATURES

- 3D gyrating motion
- 12" square aluminum tray with flat mat (standard)
- 8" square or 12" square double stacking tray (optional)
- Dimpled rubber mat (optional)
- Supplied with flat mat
- CSA approved for US and Canada
- CE marked
- 2 year warranty

ACCESSORIES

Catalog No.	Description
260288	Optional 8 inch Tray (Includes Platform, 1 each, flat and dimpled mat)
C 2202208	Dimpled rubber mat 12"
C 2202162	Flat rubber mat 12"
C 2202774	Dimpled rubber mat 8"
C 2202775	Flat rubber mat 8"
260312	Stacking tray 12"

Tube Rotator	
Catalog No.	260650
Speed	4 - 80 rpm
Carousel angle	90 - 180 degrees
Timer (with auto shut off)	Up to 99 hrs and 59 minutes
Carousel capacity	46 9 - 11 mm tubes, 28 14 - 18 mm tubes, 14 27 - 34 mm tubes, 63 18.5 mm culture tubes (drum rotor) and 63 14 mm culture tubes (drum rotator)
Product Dimensions W x D x H	304 mm x 431 mm x 381 mm (12" x 17" x 15")
Power	19.4 V dc, 474 A via 100/ 240 V AC adapter
Shipping Weight	7.25 kgs. (16 lbs.)

Boekel Orbital Shaker	
Catalog No.	270200
Speed	50 to 300 rpm
Orbit	18 mm
Timer (with auto shut off)	User defined up to 99 hours
Time setting resolution	1 min
Product Dimensions W x D x H	267 mm x 280 mm x 100 mm (10.5" x 11" x 4")
Tray Dimensions W x D	254 mm x 203 mm (10" x 8")
Power	12 V dc, 0.5 A via 115/ 230 V AC adapter
Shipping Weight	5.8 kgs. (13 lbs.)



270290



270291



270292



270293

Catalog No.	Orbitron III (115 VAC)	Orbitron III (230 VAC)
	260301	260301-2
Warranty	2 years	2 years
Operating Temperature	Can operate at 4°C to 40°C	Can operate at 4°C to 40°C
Variable speed	6-125 orbits/min	6-125 orbits/min
Variable tilt angle	0° - 23°	0° - 23°
Load limit	9 kgs. (20 lbs.)	9 kgs. (20 lbs.)
Timer	Adjustable up to 12 hours	Adjustable up to 12 hours
Product Dimensions W x D x H	30.5 cm x 33.0 cm x 22.9 cm (12" x 13" x 9")	30.5 cm x 33.0 cm x 22.9 cm (12" x 13" x 9")
Power	115 VAC; 50/60Hz, 6 watts, 0.6 amps	230VAC; 50/60Hz, 60 Watts, 0.7 amps CE marked
Shipping Weight	9.9 kgs. (22 lbs.)	9.9 kgs. (22 lbs.)