

Boekel Groovin' Tubes Compact Benchtop Shaker/Incubator

THREE INSTRUMENTS IN ONE

- Microtube thermoshaker
- Compact benchtop incubator, without shaking
- Microtube shaker without temperature control, also suitable for cold room operation

FEATURES

- Rapid heat up
- Continuous or timed operation
- Programmable timer
- 2 year warranty
- Power pac is UL, CSA and CE approved



SPECIFICATIONS

Catalog No. 270500	Groovin' Tubes - Shaker Incubator
Speed	250 to 1400 rpm
Orbit	2mm
Uniformity	±0.2°C
Uniformity	±0.1°C
Timer (with auto shut off)	1 min to 96 hrs
Time setting resolution	1 min
Temperature range	Ambient +5 to 100°C
Product Dimensions W x D x H	205 mm x 230 mm x 130 mm (8" x 9" x 5")
Power	12 V dc, 0.5 A via 115/ 230 V AC adapter

Block not included.
See interchangeable blocks listed below.

ACCESSORIES

Catalog No.	Interchangeable Platforms
270590	20 x 1.5/2.0mL microtubes
270591	20 x 0.5 mL and 12 x 1.5 mL microtubes

