


## Boekel Microjive Microplate Shakers

### FEATURES

- Compact with low profile and small footprint
- Reliable, regulated shaking for two or four microplates
- Simple screw attachments allow quick and easy fitting of any standard depth 96- or 384-well plates, holding the plates securely in place
- 2 year warranty
- Power pac is UL, CSA and CE approved

### ACCESSORIES

Catalog No.	Interchangeable Platforms
270390	Platform for four microplates



### SPECIFICATIONS

	Microjive-2 - 2 Plate Capacity	Microjive-4 - 4 Plate Capacity
Catalog No.	270300	270340
Speed	150 to 1500 rpm	150 to 1500 rpm
Orbit	2 mm	2 mm
Timer (with auto shut off)	1 min to 4 hrs	1 min to 4 hrs
Time setting resolution	10 sec	10 sec
Product Dimensions W x D x H	175 mm x 185 mm x 90 mm (7" x 7.25" x 3.5")	175 mm x 185 mm x 90 mm (7" x 7.25" x 3.5")
Power	12 V dc, 0.5 A via 115/ 230 V AC adapter	12 V dc, 0.5 A via 115/ 230 V AC adapter

## Boekel Jitterbug-2/Jitterbug-4 Microplate Thermoshaker

### FEATURES

- Accommodates two standard microplates or strip trays
- Independently set temperature, shaking speed and mixing duration
- Chamber in a chamber design for excellent stability
- Smooth, quiet operation
- Heated lid to help eliminate condensation and improve uniformity
- 2 year warranty
- Accommodates Std. 96, 384, and 536 Plates
- Power pac is UL, CSA and CE approved



### SPECIFICATIONS

	Jitterbug-2 - 2 Plate Capacity	Jitterbug-4 - 4 Plate Capacity
Catalog No.	270400	270440
Speed	250 to 1200 rpm	250 to 1200 rpm
Orbit	2 mm	2 mm
Stability	±0.1°C	±0.1°C
Uniformity	±0.2°C	±0.2°C
Timer (with Auto-Switch off)	1 min to 96 hrs	1 min to 96 hrs
Time setting resolution	1 min	1 min
Temperature range	Ambient +5 to 60°C	Ambient +5 to 60°C
Product Dimensions W x D x H	265 mm x 270 mm x 145 mm (10.5" x 10.6" x 6")	385 mm x 270 mm x 140 mm (15.5" x 15.5" x 6")
Power	12 V dc, 4.5 A via 115/ 230 V AC adapter	12 V dc, 4.5 A via 115/ 230 V AC adapter