

Boekel Microplate Shakers

FEATURES

- Compact with low profile and small footprint
- Reliable, regulated shaking for two or four microplates
- Simple screw attachments allow quick and easy fitting of any standard depth 96- or 384-well plates, holding the plates securely in place
- LCD Display
- Power pac is UL, CSA and CE approved
- CE marked
- 2 year warranty



ACCESSORIES

Catalog No.	Interchangeable Platforms
270390	Platform for four microplates

	2 Plate Capacity	4 Plate Capacity
Catalog No.	270300	270340
Speed	250 to 1000 rpm	250 to 1000 rpm
Orbit	2 mm	2 mm
Timer (with auto shut off)	1 min to 96 hrs	1 min to 96 hrs
Time setting resolution	10 sec	10 sec
Product Dimensions W x D x H	175 mm x 185 mm x 90 mm (7" x 7.25" x 3.5")	175 mm x 185 mm x 90 mm (7" x 7.25" x 3.5")
Power	12 V dc, 0.5 A via 115/ 230 V AC adapter	12 V dc, 0.5 A via 115/ 230 V AC adapter
Shipping Weight	4.5 kgs. (10 lbs.)	4.5 kgs. (10 lbs.)

Boekel Orbital Shaker

FEATURES


- Compact with low profile and small footprint
- Smooth horizontal orbital motion for mixing in bottles, flasks, and beakers
- Soft start and noise-free operation
- Can be used in an incubator 4°C to 40°C
- Platform not supplied. See interchangeable platforms listed below.
- Power pac is UL, CSA and CE approved
- CE marked
- 2 year warranty

ACCESSORIES


Catalog No.	Interchangeable Platforms
270290	6 x 250 mL flasks/400 mL beakers
270291	12 x 100 mL flasks/150 mL beakers
270292	Universal platform with adjustable bars
270293	Flat platform with rubber mat

Orbital Shaker


Catalog No.	270200
Speed	50 to 300 rpm
Orbit	18 mm
Load Limit	2.5 kgs. (5.5 lbs.)
Timer (with auto shut off)	User defined up to 99 hours
Time setting resolution	1 min
Tray Dimensions W x D	254 mm x 203 mm (10" x 8")
Product Dimensions W x D x H	267 mm x 280 mm x 100 mm (10.5" x 11" x 4")
Power	12 V dc, 0.5 A via 115/ 230 V AC adapter
Shipping Weight	5.8 kgs. (13 lbs.)




270290



270291



270292



270293