


Boekel Microplate Replicators

FEATURES

- Pins are made of stainless steel for durability, autoclaving and flame sterilization
- Economical substitute for expensive robotics systems
- Models 140500 and 140501 include an alcohol reservoir and cleaning pads for debris removal



SPECIFICATIONS

Catalog No.	Product Dimensions W x D x H	Net Weight	Shipping Weight	Pin Length	Pin Size	Pin Design
96 Pin Replicator Model 140500	11 cm x 7.6 cm x 11 cm (4.5" x 3" x 4.5")	0.6 kg (1.25 lbs)	1.5 kg (1.75 lbs)	2.54 cm (1")	1.55 mm x 1.55 mm (0.0625" x 0.0625")	Flat, square
48 Pin Replicator Model 140501	11 cm x 7.6 cm x 11 cm (4.5" x 3" x 4.5")	0.45 kg (1 lbs)	0.7 kg (1.75 lbs)	2.54 cm (1")	1.55 mm x 1.55 mm (0.0625" x 0.0625")	Flat, square
384 Pin Replicator Model 140384	12.8 cm x 8.5 cm x 12.3 cm (5" x 3.75" x 4.75")	0.6 kg (1.25 lbs)	0.9 kg (2 lbs)	2.2 cm (1")	0.019 cm dia. (0.05")	Flat, square
Template Model 140384T1	15.1 cm x 11 cm x 6.6 cm (6" x 4.25" x 2.5")	0.6kg (1.25 lbs)	0.9 kg (2 lbs)	NA	NA	NA

Boekel Tap Dancer Mini Tube Vortex Mixer

FEATURES

- Compact, with low profile and small footprint
- Pressure sensitive cup
- Continuous or touch operation
- Rubber suction pads prevent walking, absorb vibration
- Low voltage supply for safe cold room use
- 2 year warranty
- Power pac is UL, CSA and CE approved



SPECIFICATIONS

Catalog No. 270100	Tap Dancer - Vortex Mixer
Speed	300 to 6000 rpm/ continuous or burst
Orbit	2 mm
Tube diameter (maximum)	20 mm
Ambient Operating Temperature Range	+4°C to +45°C
Product Dimensions W x D x H	105 mm x 150 mm x 85 mm (4.13" x 5.9" x 3.35)
Power	12 V dc, 0.5 A via 115/ 230 V AC adapter