



### Mini Orbitron™ Platform Rotator

#### FEATURES

- Provides gentle, continuous mixing of specimens
- Extremely quiet
- Auto reverse
- Low profile and small footprint minimizes bench space
- Accommodates a variety of blood tubes, Centrifuge tubes, glass vials, or plates
- Supplied with two white, non-slip mats, one dimpled and one flat. Rubber mats can be cleaned with soap and water and autoclaved too.
- 2 year warranty

### Orbitron I™ Platform Rotator

#### FEATURES

- Gentle, low foaming agitation
- Small footprint
- Can be used in incubator or cold room
- Accepts many different sized tubes simultaneously with no need for unwieldy racks or clamps
- Thin power cord easily passes through incubator door gaskets
- Rugged steel platform has notches for bungee cords if necessary
- Fixed speed and position platform (approximately 20 orbits per minute @ 23° angle)
- Supplied with 12" x 12" tray and dimpled mat
- Power pac is UL, CSA and CE approved
- CE marked
- 2 year warranty



### Orbitron II™ Platform Rotator

#### FEATURES

- Gentle, low foaming agitation
- Small footprint
- Can be used in cold room or incubator
- Accepts a wide variety of tubes, plates, dishes, simultaneously with no need for unwieldy racks or clamps
- Thin power cord easily passes through incubator door gaskets
- Rugged steel platform has notches for bungee cords if necessary
- Variable speed of 1-23 orbits per minute
- Supplied with two white, non-slip mats, one dimpled and one flat. Rubber mats can be cleaned with soap and water and autoclaved too.
- Power pac is UL, CSA and CE approved
- CE marked
- 2 year warranty

### Boekel Rocker II™

#### FEATURES

- Ambient temperature range 4°C to 40°C
- Low voltage power is supplied by wall mount transformer
- Accepts a wide variety of flasks and plates and tubes
- Adjust tilt angle and speed without the need for tools
- Thin power cord easily passes through incubator door gaskets
- Optional stacking tray and dimpled mat (for test tubes) also available
- Supplied with flat non-skid mat
- Power pac is UL, CSA and CE approved
- CE marked
- 2 year warranty

#### ACCESSORIES

Catalog No.	Description
260355	Stacking Tray (Supplied with rubber non-skid mat)
C2202208	Dimpled Mat



	Mini Orbitron (115 VAC)	Mini Orbitron (230 VAC)
<b>Catalog No.</b>	<b>260100</b>	<b>260100-2</b>
Warranty	2 years	2 years
Fixed Speed	24 rpm	24 rpm
Fixed Tilt Angle	20°	20°
Load Capacity	0.8 kg (1.8 lbs)	0.8 kg (1.8 lbs)
Operating Temperature	+4°C to 45°C	+4°C to 45°C
Product Dimensions W x D x H	19 cm x 12.7 cm x 13.97 cm (7.5" x 5" x 5.5")	
Tray Size W x D	24.3 cm x 13.1 cm (9.75" x 5.25")	24.3 cm x 13.1 cm (9.75" x 5.25")
Net Weight	1 kg (2.21 lbs)	1 kg (2.21 lbs)
Power	115 VAC, 60 Hz, 0.2 amps	230 VAC, 50 Hz, 0.2 amps
Shipping Weight	1.3 kgs. (3 lbs.)	1.3 kgs. (3 lbs.)

	Orbitron I (115 VAC)	Orbitron I (230 VAC)
<b>Catalog No.</b>	<b>260200</b>	<b>260200-2</b>
Warranty	2 years	2 years
Temperature	4°C to 65°C (in a cold room or incubator)	4°C to 65°C (in a cold room or incubator)
Fixed Speed	Approx. 20 orbits/min	Approx. 20 orbits/min
Fixed Tilt Angle	23°	23°
Load Limit	1 kg (2.25 lbs)	1 kg (2.25 lbs)
Product Dimensions w/platform W x D x H	30.5 cm x 30.5 cm x 25.4 cm (12" x 12" x 10")	30.5 cm x 30.5 cm x 25.4 cm (12" x 12" x 10")
Product Dimensions	30.5 cm x 30.5 cm (12" x 12")	30.5 cm x 30.5 cm (12" x 12")
Shipping Weight	4.4 kg (10 lbs)	4.4 kg (10 lbs)
Power	115 VAC; 50/60Hz, 6 watts, 0.5 amps	230VAC; 50/60Hz, 6 Watts, 0.5 amps CE marked

	Orbitron II (115 VAC)	Orbitron II (230 VAC)
<b>Catalog No.</b>	<b>260250</b>	<b>260250-2</b>
Warranty	2 years	2 years
Temperature	Can operate at 4°C to 40°C	Can operate at 4°C to 40°C
Variable Speed	6-23 orbits/min	6-23 orbits/min
Variable Tilt Angle	0° - 23°	0° - 23°
Load Limit	1 kg (2.25 lbs)	1 kg (2.25 lbs)
Tray Dimensions W x D	30.5 cm x 30.5 cm (12" x 12")	30.5 cm x 30.5 cm (12" x 12")
Product Dimensions W x D x H	30.5 cm x 30.5 cm x 25.4 cm (12" x 12" x 10")	30.5 cm x 30.5 cm x 25.4 cm (12" x 12" x 10")
Shipping Weight	5.5 kg (12 lbs)	5.5 kg (12 lbs)
Power	115 VAC; 50/60Hz, 6 Watts, 0.5 amps	230VAC; 50/60Hz, 6 Watts, 0.5 amps CE marked

	Rocker II (115 VAC)	Rocker II (230 VAC)
<b>Catalog No.</b>	<b>260350</b>	<b>260350-2</b>
Tilt Angle	0° to 20°	0° to 20°
Adjustable Speed	3 - 25 rpm	3 - 25 rpm
Product Dimensions w/platform W x D x H	30.5 cm x 30.5 cm x 17.8 cm (12" x 12" x 7")	
Product Dimensions w/o platform W x D x H	19.1 cm x 17.8 cm x 12.7 cm (7.5" x 7" x 5")	
Net Weight	3.9 kg (8.5 lbs)	3.9 kg (8.5 lbs)
Shipping Weight	5.9 kg (13 lbs)	5.9 kg (13 lbs)
Stacking Tray Clearance	6.3 cm (2.5")	6.3 cm (2.5")
Stacking Tray Weight	0.7 kg (1.5 lbs)	0.7 kg (1.5 lbs)
Power	12 V dc, 0.5 A via 115/ 230 V AC adapter CE marked	
Recommended Load Limit	2.7 kg (6 lbs)	2.7 kg (6 lbs)