



Boekel Jitterbug-2/Jitterbug-4 Microplate Thermoshaker

FEATURES

- Accommodates two or four standard microplates or strip trays
- Independently set temperature, shaking speed and mixing duration
- Chamber in a chamber design for excellent stability
- Smooth, quiet operation
- Temperature calibration
- Heated lid to provide greater uniformity
- Accommodates Std. 96, 384, and 536 Plates
- Power pac is UL, CSA and CE approved
- CE marked
- 2 year warranty

	Jitterbug-2 - 2 Plate Capacity	Jitterbug-4 - 4 Plate Capacity
Catalog No.	270400	270440
Speed	250 to 1200 rpm	250 to 1200 rpm
Orbit	2 mm	2 mm
Stability	±0.1°C	±0.1°C
Uniformity	±0.2°C	±0.2°C
Timer (with auto shut off)	1 min to 96 hrs	1 min to 96 hrs
Time setting resolution	1 min	1 min
Temperature range	Ambient +5 to 60°C	Ambient +5 to 60°C
Product Dimensions W x D x H	265 mm x 270 mm x 145 mm (10.5" x 10.6" x 6")	385 mm x 270 mm x 140 mm (15.5" x 15.5" x 6")
Power	12 V dc, 4.5 A via 115/ 230 V AC adapter	12 V dc, 5.0 A via 115/ 230 V AC adapter
Shipping Weight	7.25 kgs. (16 lbs.)	7.25 kgs. (16 lbs.)

Boekel Mini Tube Rocker

FEATURES

- Provides gentle rocking action for mixing of specimens
- Low profile and small footprint minimizes bench space
- Fixed speed - 18 rpm
- Fixed tilt angle - 35°
- Maximum capacity - 8 tube tray
- Supplied with non-slip autoclavable tube pad (silicon rubber)
- Power pac is CSA, CE and PSE approved
- 2 year warranty



Catalog No. 260450	Tube Rocker
Speed	18 RPM's fixed speed
Tilt Angle	35° fixed angle
Tube Capacity	Up to 8 tubes maximum 120 mm long
Load Capacity	.8 kg (1.8 lbs.)
Operating Temperature	+4°C to +40°C
Product Dimensions W x D x H	218 mm x 120 mm x 127 mm (8.6" x 4.7" x 5")
Power	12 V dc, 0.5 A via 100/ 240 V AC adapter
Shipping Weight	2.26 kg (5 lbs.)

Boekel Mini Tube Vortex Mixer

FEATURES

- Compact, with low profile and small footprint
- Pressure sensitive cup
- Continuous or touch operation
- Rubber suction pads prevent walking, absorb vibration
- Low voltage supply for safe cold room use
- Power pac is UL, CSA and CE approved
- CE marked
- 2 year warranty



Catalog No. 270100	Vortex Mixer
Speed	300 to 6000 rpm/ continuous or burst
Orbit	2 mm
Tube Diameter (maximum)	20 mm
Operating Temperature	+4°C to +40°C
Product Dimensions W x D x H	105 mm x 150 mm x 85 mm (4.13" x 5.9" x 3.35")
Power	12 V dc, 0.5 A via 115/ 230 V AC adapter
Shipping Weight	1.8 kgs. (4 lbs.)