



## Variable Speed Nutator

### FEATURES

- Variable speed from 4-24 rpm's
- Low profile and small footprint minimizes bench space
- Quiet operation
- Provides gentle, continuous mixing of specimens
- Accommodates a variety of blood tubes and centrifuge tubes
- Supplied with one rubber dimple pad. Can be washed with mild soap and water.
- 2 Year Warranty

	Variable Speed Nutator (115 VAC)	Variable Speed Nutator (230 VAC)
<b>Catalog No.</b>	<b>201100</b>	<b>201100-2</b>
Speed	Variable 4-24 rpm	Variable 4-24 rpm
Tilt Angle   Fixed	20°	20°
Load Capacity	0.8 kg (1.8 lb)	0.8 kg (1.8 lb)
Operating Temperature	+4°C to +40°C	+4°C to +40°C
Tray Size	24.3 cm (9 3/4") W x 13.1 cm (5 1/4")D	24.3 cm (9 3/4") W x 13.1 cm (5 1/4")D
Foot Print W x D x H	18.7 cm (7-1/2") x 12.5 cm (5") x 13.7 (5-1/2")	18.7 cm (7-1/2") x 12.5 cm (5") x 13.7 (5-1/2")



## Boekel Tube Rotator

### FEATURES

- Digital microprocessor control
- Adjustable speed provides gentle or vigorous mixing action
- Adjustable angled head offers mixing angles from 90 to 180 degrees
- Five carousel options provide greater versatility for the most tube mixing needs
- Timer features offers time-down or time-up counting feature
- Includes 260653 - 14 Place Tube Carousel for 50mL
- Power pac is UL, CSA and CE approved
- CE marked
- 2 year warranty

### ACCESSORIES

Catalog No.	Interchangeable Carousel
260651	46 place tube carousel for 1.5/ 2.0 ml tubes
260652	28 place tube carousel for 15 mL tubes
260653	14 place tube carousel for 50 mL tubes
260654	63 place tube carousel for 18.5 mm cul. tubes
260655	63 place tube carousel for 14 mm cul. tubes

### Tube Rotator

Catalog No.	260650
Speed	4 - 80 rpm
Carousel angle	90 - 180 degrees
Timer (with auto shut off)	Up to 99 hrs and 59 minutes
Carousel capacity	46 9 - 11 mm tubes, 28 14 - 18 mm tubes, 14 27 - 34 mm tubes, 63 18.5 mm culture tubes (drum rotor) and 63 14 mm culture tubes (drum rotor)
Operating Temperature	+4°C to +40°C
Product Dimensions W x D x H	304 mm x 431 mm x 381 mm (12" x 17" x 15")
Power	24.0 V dc, 2.70 A via 100 - 240 V AC adapter
Shipping Weight	7.25 kgs. (16 lbs.)