



Small & Large Desiccators



Applications

- Storing reagents, IR cells or soil samples.
- Drying biological specimens
- Cooling crucibles
- Storage sludge samples

Features

- Durable, hand crafted construction. Fabricated from stainless steel and glass.
- Superior seal integrity
- Rigid door full length door hinges
- One piece neoprene door gasket
- Shelves are adjustable in half-inch increments
- Two different sizes to select from
- Automatic desiccant regeneration available

Specifications

Model	Model 1340	Model 1342	Model 1344
Overall W x L x H	31.8cm x 31.1cm x 31.1cm (12 1/2" x 12 1/4" x 12 1/2")	31.8cm x 31.1cm x 31.1cm (12 1/2" x 12 1/4" x 12 1/2")	41.9cm x 58.4cm x 57.2cm (16 1/2" x 23" x 22 1/4")
Chamber W x L x H	29.2cm x 25.4cm x 30.5cm (11 1/2" x 10" x 12")	29.2cm x 25.4cm x 30.5cm (11 1/2" x 10" x 12")	40cm x 50.2cm x 55.2cm (15 3/4" x 19 3/4" x 21 3/4")
Net Weight	10.4kg (23 lbs)	10.4kg (23 lbs)	25.9kg (57 lbs)
Shelf Info	24.8 x 20.3cm 9.75" x 8"	24.8 x 22.9cm 9.75" x 9"	48.3 x 34.3cm 19" x 13.5"

Ordering Information

Description

Small Desiccator, two cement board shelves & desiccator tray
 Small Desiccator, two stainless steel shelves & desiccator tray
 Large Desiccator, two stainless steel shelves & desiccator tray

Catalog/Model Number

1340
 1342
 1344

WARNING NOTE: No vacuum or pressure should be applied to 1340, 1342 or 1344 desiccators

Accessories

Shelves' Catalog Number

D1004369
 D1004371
 C000227

Use with

1342
 1344
 1340