



Contrast Media Warmer

Part Numbers: 133001-1 (115 volt/50-60 Hz)
133001-2 (230 volt/50-60 Hz*)

The Boekel Scientific Contrast Media Warmer with a locking door is a safe option for hospital and lab settings where security is required. Analog controls with an LCD temperature display make this 0.8-cubic-foot incubator an economical choice.

This contrast warmer comes with three removeable metal shelves and can accommodate up to four, which are available for separate purchase. The unit has a pilot light for heating indication and an external port for a thermometer or secondary temperature monitoring.



Model Numbers	133000-1 for 115V and 133000-2 for 230V
Chamber Volume	0.8 cu. ft.
Chamber Dimen. (in.)	11-1/2 x 12 x 10-1/2
W x D x H (cm)	29.2 x 30.5 x 26.7
Temperature Range	Ambient +5°C to 60°C
Stability @ 37°C	±1.5°C
Uniformity @ 37°C	±1°C
Shelves Included	2
Max. No. Shelves*	4
Shelf Area	1.7 ft.2
Dimensions (in.)	13 x 14 x 16.25
W x D x H (cm)	33 x 35.6 x 41.3
Shipping Weight	22 lb./10 kg

Features

- Locking Door for Security
- LCD temperature display
- Analog Dial Control
- Heating pilot light indicates when heater is on
- Top port for thermometer or secondary temperature probe
- Small footprint
- Magnetic latch
- 115 and 230 Volt models
- Includes 2 shelves

Contrast Media Warmer

Applications

- Pharmacy use
- Warm packs for phlebotomy
- Drinking water testing
- Water quality testing
- Incubating assays
- Throat cultures
- Incubating E. coli
- Broth incubation
- Incubating petrifilm
- Heating histology slides
- Food sample testing
- Incubating petri dishes

