





# Boekel Scientific™ Thermal Blocks

For Use with Boekel Scientific™ Heat/Cool & Heat Only Thermal Mixers II


**270810**  
0.5mL Tube Block



**270811**  
1.5mL Tube Block



**270812**  
2.0mL Tube Block




Block	Block, 24 x 0.5mL Tubes
Part No.	270810
Max Temperature (°C)	100
Max Speed (RPM)	2000
Min Speed (RPM)	0


Block	Block, 24 x 1.5mL Tubes
Part No.	270811
Max Temperature (°C)	100
Max Speed (RPM)	2000
Min Speed (RPM)	0

Block	Block, 24 x 2.0mL Tubes
Part No.	270812
Max Temperature (°C)	100
Max Speed (RPM)	2000
Min Speed (RPM)	0


**270813**  
5.0mL Tube Block



**270814**  
12mm Tube Block



**270815**  
Cryo Tube Block



Block	Block, 8 x 5.0mL Tubes
Part No.	270813
Max Temperature (°C)	100
Max Speed (RPM)	1000
Min Speed (RPM)	0

Block	Block, 24 x 12mm Tubes
Part No.	270814
Max Temperature (°C)	100
Max Speed (RPM)	2000
Min Speed (RPM)	0




Block	Block, 24 x Cryo Tubes
Part No.	270815
Max Temperature (°C)	100
Max Speed (RPM)	2000
Min Speed (RPM)	0






# Boekel Scientific™ Thermal Blocks

For Use with Boekel Scientific™ Heat/Cool & Heat Only Thermal Mixers II



**270816**  
**15mL Tube Block**



**270817**  
**50mL Tube Block**



**270818**  
**Microplate Block**






Block	Block, 8 x 15mL Conical Tubes
Part No.	270816
Max Temperature (°C)	100
Max Speed (RPM)	1000
Min Speed (RPM)	0

Block	Block, 4 x 50mL Conical Tubes
Part No.	270817
Max Temperature (°C)	100
Max Speed (RPM)	1000
Min Speed (RPM)	0



Block	Block, Standard and Deepwell Microplate	
Part No.	270818	
Max Temperature (°C)	100	
Max Speed (RPM)	Standard Deepwell	3000 1000
Min Speed (RPM)		0

**270819**  
**96 Well PCR Plate Block**



Block	Block, 96 well PCR plate
Part No.	270819
Max Temperature (°C)	100
Max Speed (RPM)	2000
Min Speed (RPM)	0

**270820**  
**384 Well PCR Plate Block**



Block	Block, 384 well PCR plate block
Part No.	270820
Max Temperature (°C)	100
Max Speed (RPM)	3000
Min Speed (RPM)	0