

HISTOLOGY SAMPLE SOLUTIONS


CASSETTE STORAGE FOR FFPE SAMPLES



Boekel Scientific

www.boekelsci.com

Boekel Scientific introduces a robust, modular storage solution for formalin fixed paraffin embedded tissue samples.



Stackable Design

Stack up to three units high with the non-tip design

(Two Units Pictured)

Durable Construction

Easy to remove crystal clear trays with metal outer housing



High Capacity

Store up to 1,200 cassettes with foam organizers included

Interchangeable Labels

Easy-slide drawers with customizable labels



Cassette Storage Unit Part Number	143000
Additional Tray Inserts Part Number	143100 (Case of 20 inserts)
Outer Dimensions	W 13.25" x D 13.5" x H 7.75"
Cassette Capacity	Approximately 1,200 Cassettes per Unit
Unit Stacking Capacity	3 Shelving Units High
Shelf Material of Construction	Polycarbonate
Housing Material of Construction	Epoxy Painted Metal



855 Pennsylvania Boulevard | Feasterville, PA 19053 USA | 800.336.6929

p. 215.396.8200 | f: 215.396.8264 | e-mail: boekel-info@boekelsci.com

©Copyright, Boekel Scientific, 2013